

Snap-Fit™

GAC Molar Bands



Prewelded Molar Bands

Various attachments available

Snap-Fit Molar Bands

Our Snap-Fit Molar Bands offer greater elasticity and superior fit. GAC's exclusive process produces bands with improved morphology, elasticity and memory. Snap-Fit bands are created with high contouring and snap right into place when the right size is implemented. You simply press on the disto-lingual and then the mesio-labial and the band "snaps" into place. After fitting and removal, the band returns to its original size, assuring that the cemented band is as accurate as the initial fit.

Snap-Fit Molar Bands are offered in graduated, proportional sizes, maxillary and mandibular, left and right (32 sizes for first molars, 30 sizes for second molars). Snap-Fit Molar Bands are available plain or pre-welded.



Snap-Fit are the choice of many leading orthodontists world wide. A quarter of a century of experience has proven that they can dramatically reduce your loose band rate if properly sized. It is no wonder that

Snap-Fit bands have been well accepted by thousands of clinicians across the world.

All Snap-Fit Bands are marked using our patented Electro Chemical Etching process for easy-to-read permanent identification, and interiors are micro-etched for roughness to increase bond strength.

Micro etched interiors
to increase bond strength

Polished exteriors
for maximum patient comfort



Guaranteed weld strength
for extra peace of mind

Electro-chemically etched
sizes for easy ID